



MSAB Tablet

The MSAB Tablet is a lightweight and easy to use turnkey solution. It is designed for the frontline, where mobile data recovery needs to happen quickly to ensure timely access to evidence.

With a touch screen interface, the Tablet is designed to quickly and easily recover data from mobile devices. It allows you to connect the mobile device, touch the screen and extract the data within minutes.



Fast extraction times



User Defined Workflow and Case Fields



Simplified user interface



Portable and flexible

Rapid Access to Mobile Data

With its ease of use the Tablet provides first responders with real-time evidence and intelligence gathering capability in the most flexible form factor to date. Whilst complex cases or devices will always need highly trained personnel, the Tablet can help speed up the 'triage' process of assessment by allowing investigators quick access.

Extract and view

The Tablet offers support for extraction of mobile devices through cable and Bluetooth. Users perform extractions in a simplified extraction wizard, in an easy to use touch based interface. Extracted data can be overviewed and explored in the Viewer application. Through a powerful search function, and easy to use graphics, specific selections of the extracted data can be displayed. The extraction creates a secure XRY file, containing all extracted data from the device.

Workflow Configuration

By the creation of workflow sequences a user can be guided through a tailored extraction process. The workflow can be configured for your specific organizational needs to include or exclude tasks that meet your policies. This simplifies and secures the process of device extraction and examination. A typical workflow sequence can also include help screens to guide the user in the extraction process.

User Administration

User management offers control of tasks being performed on the Tablet. When a user is created the administrator assigns rights to workflow sequences for that user. Examples of workflow sequences include Extract data, View extraction, SIM id cloner as well as custom made multistep workflow sequences. In standard user mode the device is locked down and all user activities are logged. This way activities are traceable and the extraction process is transparent. User Administration functions also include privileges to delete and export all user files.

Kit contents

- XRY Logical and/or Physical 12 month License with updates
- XAMN Viewer software for analysis and reporting
- 11th Gen Intel Core I7 Processor tablet with anti-glare and anti-smudge display
- Ruggedized Case + detachable keyboard
- XRY Cable Kit
- AC/DC Adapter
- Memory Card Reader - Write Protected
- SIM Card Reader
- Free cable updates for 12 months

Specifications

- 11th Gen Intel Core I7 Processor
- 13" Touch Screen
- 16GB RAM
- 1TB Solid State Drive
- Win 10 64Bit OS
- 3 Year warranty
- Tablet Unit: 288.4 x 207.9 x 8.44 mm (11.35" x 8.18" x 0.33")
- Keyboard: 284.2 x 215.6 x 5.1 mm (11.18" x 8.48" x 0.2")
- Weight: 800g (1.8 lbs) (Tablet only)



Contact sales@msab.com to get started!

MSAB

Trusted Partner in Digital Forensics